

# Multifunctional Artificial Teeth & Toothbrush Cleaner

Thank you for purchasing our “Multifunctional Artificial Teeth & Toothbrush Cleaner”. Compared with the other current denture cleaners, our product includes UV sterilization light, nano-silver technology with Ultrasonic wave, and high speed vibration to easily get rid of debris on dentures, jewelry, toothbrush... etc. Moreover, small size is perfect for home and traveling. By carefully following the below instruction, you will be able to have trouble-free operation to enjoy great oral hygiene for whole life.

## **Function**

1. Our product is applying the newest vibration and sterilization effect in cleaning dentures, jewelry, toothbrush, and so on. Its vibration reaches 18,000 times per minute in order to power away stains from dentures, toothbrush, and jewelry. In addition, we adopt the latest Japanese techniques of UV light sterilization and nano-silver mixture container to deodorize, sterilize and wash at one step.
2. This multi-function washer combines with the latest electronic sterilization technology and nano-silver bacteriostasis. Specially, no consumable replacement is needed.
3. Build-in power save timer shuts the power automatically after 20 minutes when you turn it on.
4. You do not have to use any other chemical detergent but water. Through the latest UV light sterilization and nano-silver bacteriostasis ultrasonic wave-hygienic, it can reach 100% spotless effect without residual chemical which might damage your denture, toothbrush, or jewelry.
5. To avoid the user see the UV light directly, the power will be off and all functions will be stopped after the main top is opened.
6. User-friendly design that the main top can magnetically close up when the users close it

## **Specification**

1. Dimension : 104 (L) x 127 (W)x 100 (D) (mm) / 4.09 (L) x 5.00 (W) x 3.94 (D) (in)
2. Weight : 360g
3. Vibration : 18,000 times per minute
4. Power Supply : AA size battery x 3 pcs or switching type AC adapter ( Universal from 100~240V )

## **Illustration**



- |   |                           |                  |                        |
|---|---------------------------|------------------|------------------------|
| 1. UV + Ultrasonic + Vibration                | 2. LED indicator          | 3. UV light only | 4. Battery Compartment |
| 5. Nano-silver powder mixture water container | 6. Toothbrush insert hole | 7. DC jack       |                        |

# Set down the denture or jewelry and then choose either one function. (Ultrasonic Washing & UV or UV only)Or just open the transparent cover on the top of the unit and put in the toothbrush as image shows.

It turns off automatically after 5 minutes when you press “ UV + Ultrasonic” to start, but UV light will stay longer for extra 15 minutes.

It turns off automatically after 20 minutes when you press “UV only” to start.

# In case you want to stop working in advance, press the function button for 3 seconds.